

Products

– Security & Sensors

SPM915

In-bed/Off-bed Monitoring



Description

In-bed/off-bed Monitoring Pad is mainly used to detect the status of the elderly or disabled patients in/out of bed. If the elderly or disabled patients leave the bed, it will notify your family and caregivers immediately through the APP.

Main Features

- ZigBee 3.0
- Low power consumption

Main Specifications

Connectivity	
ZigBee	<ul style="list-style-type: none">• 2.4GHz IEEE 802.15.4
ZigBee Profile	<ul style="list-style-type: none">• ZigBee 3.0
RF Characteristis	<ul style="list-style-type: none">• Operating frequency: 2.4GHz• Range outdoor: 100m (Open area)• Internal PCB Antenna
Physical Specifications	
Battery	<ul style="list-style-type: none">• DC 3V (Two AAA batteries)
Operating environment	<ul style="list-style-type: none">• Temperature: -10°C ~+55°C• Humidity: ≤ 85% non-condensing
Mounting Type	<ul style="list-style-type: none">• Lay the pressure sensing pad horizontally on the bed